

# ANYLAB F1

## High-Sensitive IVD POCT System

ANYLAB-F1 is a POCT (Point Of Care Test) system that is based on immunofluorescence and test diagnostic items with the samples - serum, plasma, whole blood or urine like a large centralized system, so it can obtain results easily and quickly(3~15mins). In addition, the system has used the immunofluorescence method which is 10~100 times more sensitive than the colorimetry method, so it can achieve quantitative results with a higher degree of accuracy, making the test more useful.



**5** Years Warranty

Authorized Distributor



**Corporate Office:**

Rupayan Taj 1, E-5, 5th Floor, 1/1 Box Culvert Road  
Naya Paltan, Dhaka-1000, Bangladesh, Web: [www.bio-techbd.com](http://www.bio-techbd.com)  
Phone No: 02-9352208, Mobile : +88 01711-947754  
E-mail: [suzonahamed87@gmail.com](mailto:suzonahamed87@gmail.com), [contact@bio-techbd.com](mailto:contact@bio-techbd.com)

## SPECIFICATION

Dimensions	305 x 260 x 115 mm
Weight	2.2 Kg
Power Supply	220 V, 5a
Data Presentation	LCD/Printer/Android
Working Temp.	15 ~ 35°C
Working Humidity	Max 70%
Expiration date	24months at 2-8°C (Buffer) 24months at 2-30°C (Cartridge)

Made in Korea